

Birntaster BT-B. Push button BT-B.



- Ruftaste mit Schwesternsymbol • integriertes Finde- und Beruhigungslicht • spritzwasser dichtes Gehäuse • antimikrobielle Folientastatur.

- Call button with nurse symbol • integrated finder and reassurance light • splash resistant case • antimicrobial membrane keyboard • integrated Switch 100BaseTX.

Beschreibung

Für die Bedienung durch den Patienten vom Krankbett aus, bestehend aus:

- Stirnseitig angebrachte Ruftaste mit Schwesternsymbol mit integriertem Finde- und Beruhigungslicht,
- Folientastatur mit integrierten LED's, bestehend aus:
 - 1 Ruftaste (rot),
 - 2 Lichttasten
- 2,80 Meter Anschlusskabel mit RJ45 Stecker
- Anschlusskabel mit selbstabwerfendem Stecker bei Zugbelastung.


Description

For operation by patients from their patient bed, consisting of:

- Call button with nurse symbol with integrated finder and reassurance light on the top end of the unit,
- Membrane keyboard with integrated LED's, consisting of:
 - 1 call button (red),
 - 2 lighting keys
- 2.80m connection cable with a RJ45 plug
- Connection cable with selfrelease plug which is protected against break by an unintentional unplug.

Birntaster BT-B. Push button BT-B.

Abmessungen Dimensions



Alle Maße in mm
All dimensions in mm

Technische Daten

Versorgungsspannung:	24V (20V bis 30V)
Stromverbrauch:	max: 200mA/20V (=4,0W) typ: 71mA/24V (=1,7W)
Stecker:	RJ45
Kabelzugentlastung:	max: 200N (geräteseitig)
Umgebungstemperatur:	0 bis +40°C
Relative Luftfeuchtigkeit:	5 bis 100%
Luftdruck:	≥80 kPa, bis 2000m über Meer
Schutzart:	IP32, VDE0834 Umweltklasse II
EMV:	EN50081-1 Störemission für Wohnbereich EN50082-2 Störfestigkeit für Industrie
Produktsicherheit:	EN60950
Funkentstörung:	EN55022B
Gehäuse Material:	Kunststoff ABS
Gehäuse Farbe:	weiß RAL 9016
Abmessungen:	96 x 64 x 23 mm (HxBxT)

Technical data

Supply voltage:	24V (20V up to 30V)
Current consumption:	max: 200mA/20V (=4.0W) typ: 71mA/24V (=1.7W)
Connector:	RJ45
Strain relief:	max: 200N (device-side)
Ambient temperature:	0 up to +40°C
Relative humidity:	5 up to 100%
Air pressure level:	≥80 kPa, up to 2000m above sea-level
Protection class:	IP32, VDE0834 environmental Class II
EMC:	EN50081-1 emission for residential use EN50082-2 emission for industrial environments
Product safety:	EN60950
Radio suppression:	EN55022B
Case material:	plastic ABS
Case color:	white RAL 9016
Dimensions:	96 x 64 x 23 mm (HxWxD)

Weitere Informationen sind der Technischen Dokumentation zu entnehmen.
Further informations can be obtained from the technical documentation.

SCHRACK SECONET AG

A-1120 Wien, Eibesbrunnengasse 18 • Tel.: +43-1-81157-0 • office@schrack-seconet.com

HEALTH CARE

www.schrack-seconet.com

SCHRACK
S E C O N E T

© Schrack Productions - K-DB-0028DE-EN_BT-B_V1.3
Technische Änderungen vorbehalten - subject to technical modifications